

### Technical Data Sheet

# ACS Material Analytical Testing of

# X-ray Diffraction (XRD)

(Powder/Film/Bulk)

#### **Contact Information:**

ACS Material, LLC

Address: 959 E Walnut St., Suite 100

Pasadena, CA 91106, USA

Phone: (866)-227-0656

Fax: (781)-518-0284

E-Mail: contact@acsmaterial.com

Revision: 112020

Test Service	X-ray Diffraction (XRD)
Instrument Model	X'Pert PRO MPD, D8, Ultima IV, Smartlab
Test Angle	Regular: 5~90°
	Small Angle: 0.5~10°
Test Rate	Regular: 10°/min, 5°/min, 2°/min
	Small Angle: 1°/min, 0.5°/min
Sample Type	Powder/Film/Bulk
Lead Time	2-3 weeks

#### **Notes:**

- 1. Sample Requirements:
  - Powder:Please provide 100 mg. If the amount is limited, please contact us. Please grind the sample
    well before delivering.
  - Bulk: Size within  $2 \times 2 \times 0.5$  cm. Please mark the test surface and ensure that the test surface is smooth and clean.
  - Film: Please contact us for more information.
- 2. Reference only and it may vary from sample to sample.
- 3. All samples are non-returnable and contact us if you have special requests.

**Disclaimer:** ACS Material, LLC believes that the information in this Technical Data Sheet is accurate and represents the best and most current information available to us. ACS Material makes no representations or warranties either express or implied, regarding the suitability of the material for any purpose or the accuracy of the information contained within this document. Accordingly, ACS Material will not be responsible for damages resulting from use of or reliance upon this information.